# SEEKER<sup>®</sup> Crossing Support Catheter


ENGLISH

## INSTRUCTIONS FOR USE

**CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.**

### Device Description

The SEEKER<sup>®</sup> Crossing Support Catheter (Figure 1) is a catheter consisting of an over-the-wire catheter with a fixed HUB. Radiopaque markers delineate the distal tip of the catheter and aid in catheter placement as well as lesion length measurement. The catheter includes an atraumatic tip to facilitate advancement of the catheter to and through the peripheral vasculature stenosis and chronic total occlusions. The catheter contains a hydrophilic coating that is intended to provide lubricity.

The over the wire catheter is compatible with .014", .018" and .035" guidewires and is available in 65, 90, 135, and 150 cm working lengths. The proximal portion of the catheter includes a luer hub connected to the catheter lumen. The distal portion of the catheter includes radiopaque markers that are 1 cm in length and are positioned 1 cm apart totaling 30 cm. Double markers delineate the 10cm and 20cm radiopaque markers. There is a 1mm long radiopaque marker band located 5-6mm from the distal tip. These products are not made with natural rubber latex.

Reference Table 1 for SEEKER<sup>®</sup> Product Dimensions.

Table 1

SEEKER <sup>®</sup> Product Dimensions	SK13514	SK15014	SK9018	SK13518	SK15018	SK6535M	SK9035M	SK13535M	SK15035M
Maximum guidewire, Inch	.014	.014	.018	.018	.018	.035	.035	.035	.035
Catheter working length, cm	135	150	90	135	150	65	90	135	150
Tip outer diameter, Inch	.017	.017	.021	.021	.021	.038	.038	.038	.038
Distal outer diameter, Inch	.029	.029	.033	.033	.033	.053	.053	.053	.053
Proximal outer diameter, Inch	.039	.039	.045	.045	.045	.064	.064	.064	.064
Minimum guide catheter, FR.	4	4	4	4	4	5	5	5	5
Minimum introducer sheath, FR.	4	4	4	4	4	5	5	5	5

### Indications for Use

The SEEKER<sup>®</sup> Crossing Support Catheters are recommended for guidewire exchange and infusion devices designed for use in the peripheral vascular system. The catheters are intended to support a guidewire during access into the vasculature, allow for guidewire exchanges, and provide a passage for delivery of saline solutions and/or diagnostic contrast agents.

### Contraindications

None known.

### Warnings

1. Contents supplied STERILE using ethylene oxide (EO). Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Single patient use only. Do not reuse, reprocess or resterilize.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
3. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation.
5. While retracting the catheter, if resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or catheter separation.
6. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.
7. Exceeding the maximum injection pressure of 300 psi may result in catheter failure.

### Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The SEEKER<sup>®</sup> Crossing Support Catheter shall only be used by physicians trained in the use of vascular diagnostic catheters.
3. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the crossing support catheter through a smaller size introducer sheath than indicated on the label.
4. Only use guidewires of the recommended diameter as noted on the packaging.
5. Flush catheter with sterile saline until saline is seen exiting distal tip.
6. If resistance is felt during post procedure withdrawal of the catheter, it is recommended to remove the catheter, introducer sheath and guidewire (as necessary) as a single unit.
7. Do not continue to use the catheter if the shaft has been kinked.
8. Prior to re-insertion through the introducer sheath, wipe the catheter clean with wet gauze, flush with sterile saline, and check the tip for damage and kinking.
9. Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date.
10. If the catheter is to be used as an infusion catheter, reference Table 2 to ensure proper infusion rates/pressures are utilized and not exceeded.

Table 2

SEEKER <sup>®</sup> Infusion Flow Rates (ml/sec)	Sterile Saline		Contrast*			
	150 psi	300 psi	150 psi	300 psi		
SK13514	.014	135	1.0	1.6	0.2	0.6
SK15014	.014	150	0.7	1.2	0.2	0.4
SK9018	.018	90	1.4	2.3	0.6	1.5
SK13518	.018	135	1.5	2.6	0.5	1.1
SK15018	.018	150	1.3	2.3	0.5	1.2
SK6535M	.035	65	8.2	12.0	5.9	8.8
SK9035M	.035	90	7.5	10.4	5.1	7.7
SK13535M	.035	135	6.4	8.8	4.3	6.7
SK15035M	.035	150	6.0	8.6	4.0	6.5

Flow rates based on a 75/25 contrast/saline mixture at a viscosity of 4.2 cP.

### Potential Adverse Reactions

The complications which may result from a crossing support catheter procedure include:

- Vessel dissection, perforation, rupture or total occlusion
- Embolism
- Hypotension/hypertension
- Hemorrhage or hematoma
- Infection
- Vessel Spasm
- Distal embolization of thrombotic or atherosclerotic material
- Death

### Directions for Use

#### Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources. Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date.

Do not use if packaging is damaged or opened.

#### Equipment Required

- Sterile saline solution
- Luer lock syringe
- Appropriate introducer sheath and dilator set
- Recommended guidewire

### SEEKER<sup>®</sup> Crossing Support Catheter Preparation

1. Verify the crossing support catheter size is suitable for the procedure and the selected accessories accommodate the catheter as labeled. Remove catheter from package.
2. Prepare the lumen of the catheter by attaching a syringe to the luer hub and flush with sterile saline solution.
3. Activate the hydrophilic coating by exposing the outside of the catheter to sterile water or saline.

### Use of the SEEKER<sup>®</sup> Crossing Support Catheter

1. Place the distal tip of the SEEKER<sup>®</sup> Crossing Support Catheter over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath and over the wire to the site of the lesion. If the stenosis cannot be crossed with the desired catheter, exchange guidewire for smaller diameter guidewire first, then exchange the SEEKER<sup>®</sup> Crossing Support Catheter for the smaller wire size chosen.
3. Once the wire and catheter have crossed the lesion, remove the catheter by gently pulling back on the catheter through the introducer sheath keeping your guidewire in place.

### SEEKER<sup>®</sup> Crossing Support Catheter Reinsertion

**Precaution:** Do not continue to use the SEEKER<sup>®</sup> Crossing Support Catheter if the shaft has been damaged or kinked.

**Precaution:** Prior to re-insertion through the introducer sheath, the catheter should be wiped clean with gauze and flushed with sterile saline.

1. Continue the procedure according to the "Use of the SEEKER<sup>®</sup> Crossing Support Catheter" section herein.

**Warning:** After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

### Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country. An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Copyright © 2017 C. R. Bard, Inc. All rights reserved. Printed in U.S.A.**
Copyright © 2017 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright © 2017 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in USA.
Copyright © 2017 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2017 C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Copyright © 2017 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.