# Made to fit Your Patients. Made to fit You.




## **Evolving with You**

Based on over 10 years of clinical experience, E-vita Open Neo is the next generation hybrid stent graft system for aortic arch and descending thoracic aorta repair with Frozen Elephant Trunk technique.

2005
1982

E-vita Open

Conventional Elephant Trunk 2008




Enhanced Frozen Elephant Trunk Device

E-vita Open Plus

2012



Branched FET Device



E-vita Open Neo Next generation Frozen Elephant Trunk Device

### E-vita Open Neo: Versatility Re-engineered

Supraortic trunks management with **island technique** or **sequential anastomoses** for optimized arch reconstruction

Multiple designs for **collar** anastomosis from **Zone 0 to Zone 3** depending on patient need and preferred technique

Optimized graft diameters with a wide range of stent graft options that **adapt to different anatomies and diseases** 

Redesigned delivery system for **intuitive and precise deployment** 

**Shapeable shaft or guide wire** compatibility, adjusting to your preference and needs

First Frozen Elephant

**Trunk Device** 

### Adapting to Patient's Disease

Each disease requires a unique technique and oversizing strategy, and based on this E-vita Open Neo has been created with three stent graft configurations and dedicated designs of the vascular and stent graft sections.

The three diameters of the vascular section allow for reproducible anastomosis, while the full range of options of the stent graft section provides the adequate oversizing for both aneurysms and dissections.



### STRAIGHT

Island Technique

Collar Anastomosis in Zone 2/3



#### **BRANCHED**

Sequential Anastomoses


Collar Anastomosis in Zone 1/2/3



### **TRIFURCATED**


Sequential Anastomoses

Collar Anastomosis in Zone 0/1



- \* branches for supraortic trunks available only on the Branched and Trifurcated configurations
- \*\* please refer to the IFU for the specific length of each configuration and diameter

### **Control in Your Hands**



### **Deployment Steps**



Remove protective flap or pre-shape shaft



2 Device positioning



Release trigger to unlock the system, and retract handle for progressive distal-proximal stent graft deployment



4 Complete handle retraction and withdraw the delivery system

## **Ordering Information**

#### **Straight Configuration**

Catalog Number	Vascular Graft Part Ø (mm)	Stent Graft Part Ø (mm)	Length Stent Graft (mm)
95HG2620L120-C01	26	20	120
95HG2624L120-C01	26	24	120
95HG2624L175-C01	26	24	175
95HG2626L120-C01	26	26	120
95HG2828L120-C01	28	28	120
95HG2828L180-C01	28	28	180
95HG3030L120-C01	30	30	120
95HG3030L180-C01	30	30	180
95HG3033L130-C01	30	33	130
95HG3033L180-C01	30	33	180
95HG3036L130-C01	30	36	130
95HG3036L180-C01	30	36	180
95HG3040L130-C01	30	40	130
95HG3040L180-C01	30	40	180

Sizes in italics are available on demand

#### **Branched Configuration**

Catalog Number	Vascular Graft Part Ø (mm)	Stent Graft Part Ø (mm)	Length Stent Graft (mm)
95HG2622L120-C02	26	22	120
95HG2624L120-C02	26	24	120
95HG2624L175-C02	26	24	175
95HG2626L120-C02	26	26	120
95HG2626L180-C02	26	26	180
95HG2828L120-C02	28	28	120
95HG2828L180-C02	28	28	180
95HG3030L120-C02	30	30	120
95HG3030L180-C02	30	30	180
95HG3033L130-C02	30	33	130
95HG3033L180-C02	30	33	180
95HG3036L130-C02	30	36	130
95HG3036L180-C02	30	36	180
95HG3040L130-C02	30	40	130
95HG3040L180-C02	30	40	180

#### **Trifurcated Configuration**

Catalog Number	Vascular Graft Part Ø (mm)	Stent Graft Part Ø (mm)	Length Stent Graft (mm)
95HG2624L175-C03	26	24	175
95HG2626L180-C03	26	26	180
95HG2828L180-C03	28	28	180
95HG3030L180-C03	30	30	180
95HG3033L180-C03	30	33	180
95HG3036L180-C03	30	36	180
95HG3040L180-C03	30	40	180


Vascular graft length: 100 mm without tension on every configuration and size

Perfusion branch: diameter of 10 mm and length without tension of minimum 100 mm on every configuration and size


#### **Branches Specification**

	Diameter	Length	
Branch BCT	12 mm	min. 100 mm	
Branch LCCA	8 mm		
Branch LSA	10 mm		

Branch LCCA



Branch LCCA



# **ARTIVION**

#### Learn more at artivion.com

1. M. Youssef et. al (2018) - A Multicenter Experience With a New Fenestrated-Branched Device for Endovascular Repair of Thoracoabdominal Aortic Aneurysms, 1. M. Youssel et. at (2018) - A Multicenter Experience with a New Fenestratea-Eranchea Device for Endovascular Repair of Thoracoadaominal Aortic Aneurysms, Journal of endovascular therapy, DOI: 10.1177/1526602817752147 2. A.Katsargyris et.al (2018) - Early Experience with the Use of Inner Branches in Endovascular Repair of Complex Abdominal and Thoraco-abdominal Aortic Aneurysms, European Journal of vascular and endovascular surgery, DOI: 10.11016/j.ejvs.2018.01.024 3. V. Bilman, T. Cambiaghi, A. Grandi, N. Carta, G. Melissano, R. Chiesa, L. Bertoglio (2020) - Anatomical feasibility of a new off-the-shelf inner branch stent graft (E-nside for endovascular treatment of thoraco-abdominal aneurysms, European Journal of Cardio-Thoracic Surgery, Volume 58, Isssue 6, Pages 1296-1303, https://doi.org/10.1093/ejcts/ezaa276

All products and indications are not available/approved in all markets. All trademarks are owned by Artivion, Inc. or its subsidiaries. On-X Life Technologies, Inc. Jotec GmbH and Ascyrus Medical GmbH are wholly owned subsidiaries of Artivion, Inc. © 2022 Artivion, Inc. All rights reserved.

**Artivion, Inc.** 1655 Roberts Blvd., NW Kennesaw, GA 30144 USA Phone: +1-888-427-9654 Fax: +770-590-3753

E-mail: inquiries@artivion.com For contact information by region, please visit. www.artivion.com/contact

JOTEC GmbH, Lotzenäcker 23, 72379 Hechingen, Germany