



# NAVIGATE. CAPTURE. SUCCEED.



**EN Snare<sup>®</sup>**

Endovascular Snare System

**ONE Snare<sup>®</sup>**

Endovascular Snare System

# EN Snare®

Endovascular Snare System

The only delivery catheter with **reinforced strain-resistant hub and peel-away insertion tool** on the market.

EN Snare® loop tips stay as much as **58% closer** to vessel walls than the competition.<sup>1</sup>

### Embedded Radiopaque Marker Band

identifies the snare delivery catheter's precise location under imaging

### Interlaced Nitinol Loops

designed to resist prolapse and improve capture by allowing the snare to move as a single-unit through the vessel

### Improved Vessel Navigation


facilitated by 15 degree curved tip on 6F and 7F snare delivery catheters

### Enhanced Visualization

provided by platinum strands incorporated into the loops

### High Performance Snare Delivery Catheter

designed for strength and to minimize the potential for kinking and buckling




## REINFORCED CATHETER HUB

**NEW** flexible, reinforced, strain relief hub designed to reduce buckling.

**ONE Snare®**  
2.4X more buckle force<sup>1</sup>

**EN Snare®**  
1.9X more buckle force<sup>2</sup>



# ONE Snare®

Endovascular Snare System

ONE Snare® delivery catheter provides **2.4X greater buckling resistance force** than the competition.<sup>1</sup>



**Designed for Better Foreign Body Capture** with 90 degree angle loop


**Excellent Support** from Nitinol core wire that provides shape-memory

**Kink Resistant Catheter** designed to improve strength

**Nitinol & Gold-Plated Tungsten Loop** retains snare's circular shape to support capture of foreign bodies, while allowing superior visibility.



**High Visibility** under fluoroscopy due to gold-plated tungsten loop




**UNIQUE** peel-away insertion tool for quick and efficient snare loading with EN Snare® and ONE Snare.®

**ONE Snare**  
Peel-Away insertion tool



*Competitor's non-peelable introducer*



**EN Snare**

# EN Snare®

Endovascular Snare System



## EN SNARE ORDERING INFORMATION

Catalog Number	Loop Diameter	Snare Length	Catheter Size	Catheter Length
<b>Standard Snare Kit</b>				
EN2006010	6-10 mm	120 cm	6F	100 cm
EN2006015	9-15 mm	120 cm	6F	100 cm
EN2006020	12-20 mm	120 cm	6F	100 cm
EN2007030	18-30 mm	120 cm	7F	100 cm
EN2007045	27-45 mm	120 cm	7F	100 cm
<b>Mini Snare Kit</b>				
EN1003004	2-4 mm	175 cm	3.2F	150 cm
EN1003008	4-8 mm	175 cm	3.2F	150 cm

# ONE Snare®

Endovascular Snare System



## ONE SNARE ORDERING INFORMATION

Catalog Number	Loop Diameter	Snare Length	Catheter Size	Catheter Length
<b>Standard Snare Kit</b>				
ONE500	5 mm	120 cm	4F	100 cm
ONE1000	10 mm	120 cm	4F	100 cm
ONE1500	15 mm	120 cm	6F	100 cm
ONE1501*	15 mm	80 cm	4F	65 cm
ONE2000	20 mm	120 cm	6F	100 cm
ONE2500	25 mm	120 cm	6F	100 cm
ONE3000	30 mm	120 cm	6F	100 cm
ONE3500	35 mm	120 cm	6F	100 cm
<b>Petite Snare Kit</b>				
ONE1001	10 mm	65 cm	4F	48 cm
ONE2501	25 mm	65 cm	6F	48 cm
<b>Microsnare Kit</b>				
ONE200	2 mm	175 cm	2.3-3F	150 cm
ONE201	2 mm	200 cm	2.3-3F	175 cm
ONE400	4 mm	175 cm	2.3-3F	150 cm
ONE401	4 mm	200 cm	2.3-3F	175 cm
ONE700	7 mm	175 cm	2.3-3F	150 cm
ONE701	7 mm	200 cm	2.3-3F	175 cm
<b>Snare Replacement Catheters</b>				
ONE4000	-	-	4F	100 cm
ONE6000	-	-	6F	100 cm

\*Not CE marked

Footnotes:

<sup>1</sup>Data on file, Merit 2012.

<sup>2</sup>Data on file, Merit 2014.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733